
RawCOPPER data format


Dec. 16, 2013 (svn rev.7974) Satoru Yamada

1, RawDataBlock object (to handle Raw data from COPPER board)



Source code:


https://belle2.cc.kek.jp/svn/trunk/software/rawdata/dataobjects/




Overview of RawCOPPER format (one data block from a COPPER board)

- RawCOPPER header
 - COPPER header
 - B2link HSLB header (slot A FINNESSE)
 - B2link FEE header(slot A FINNESSE)
 - Data contents(Detector buffer) (slot A FINNESSE)
 - B2link FEE trailer (slot A FINNESSE)
 - B2link HSLB trailer (slot A FINNESSE)
 - B2link HSLB header (slot B FINNESSE)
 - B2link FEE header(slot B FINNESSE)
 - Data contents(Detector buffer) (slot B FINNESSE)
 - B2link FEE trailer (slot B FINNESSE)
 - B2link HSLB trailer (slot B FINNESSE)
 - B2link HSLB header (slot C FINNESSE)
 - B2link FEE header(slot C FINNESSE)
 - Data contents(Detector buffer) (slot C FINNESSE)
 - B2link FEE trailer (slot C FINNESSE)
 - B2link HSLB trailer (slot C FINNESSE)
 - B2link HSLB header (slot D FINNESSE)
 - B2link FEE header(slot D FINNESSE)
 - Data contents(Detector buffer) (slot D FINNESSE)
 - B2link FEE trailer (slot D FINNESSE)
 - B2link HSLB trailer (slot D FINNESSE)
 - COPPER trailer
- RawCOPPER trailer

2, "RawCOPPER header" and trailer format: 2013/8/26 (Not yet confirmed)



3, COPPER header and trailer from Belle document



4, B2link FEE header/Trailer, B2link HSLB header/Trailer

Data format (Final?)

From Nakao-san's B2GM slides:

http://kds.kek.jp/getFile.py/access?contribId=143&sessionId=38&resId=0&materialId=slides&confId=13911

The format used at the telescope test

```
B2link HSLB header
HSL: 0xFFAA(16) --- B2L header | HSLB-tag(16)
B2L: '0'(1) | TT-ctime(27) | TT-type(4)
B2L: TT-tag(32)
                                                  B2link FEE header
B2L: TT-utime(32)
B2L: TT-exprun(32)
B2L: '0' | B2L-ctime(27) | debug-flag(4)
FEE: Data #0 (32)
FEE: Data #1 (32)
FEE:
FEE: Data #n (32)
                                                   B2link FEE trailer
B2L: TT-tag(16) | B2L-checksum(16)
                                                   B2link HSLB trailer
HSL: 0xFF55(16) | HSLB checksum(16)
```

tag (event number) and utime to be increased to 32-bit (done),
 HSLB-checksum, B2L-checksum to be added

5, Example: how to get information of RawCOPPER header

You can get event # info from RawCOPPER object like this;

```
StoreArray<RawCOPPER> raw_cprarray;
                                                         When StoreArray is used
for (int i = 0; i < raw_cprarray.getEntries(); i++) {</pre>
            for (int j = 0; j < raw cprarray[i].GetNumEntries(); j++) {
            Get Event number
//
                        unsigned int event_no = raw_cprarray[ i ].GetEveNo( j );
            Get RawCOPPER data block
II
                        int* buf = raw cprarray[i].GetBuffer(j);
            See contents of a data block (from RawCOPPER header to RawCOPPER trailer)
//
                       for( int k = 0; k < raw_cprarray[i].GetBlockNwords(j); k++){
                                    printf("%d\n", buf[ k ] );
            Get Detector Buffer (raw data from detector electronics)
//
                       int* buf slot a = raw cprarray[i].Get1stDetectorBuffer(i);
                       int* buf_slot_b = raw_cprarray[i].Get2ndDetectorBuffer(j);
                       int* buf_slot_c = raw_cprarray[ i ].Get3rdDetectorBuffer( j );
                       int* buf slot d = raw cprarray[i].Get4thDetectorBuffer(i);
                       int* buf_slot[4]; for( int k = 0; k < 4;k++){ buf_slot[ k ] = raw_cprarray[ i ].GetDetectorBuffer(j,k) }
            See contents of raw data from detector
                       for( int k = 0; j < raw_cprarray[ i ].Get1stDetectorNwords( j ); k++ ){
                                    printf("%d\u00e4n", buf_slot_a[k]);
                       for( int k = 0; j < raw_cprarray [ i ].Get2ndDetectorNwords( j ); k++ ){
                                    printf("%d\footnote\n", buf slot b[k]);
```

Test program to read RawCOPPER(RawCDC) data

== Detector Buffer(FINESSE A)

== Detector Buffer(FINESSE A)

== Detector Buffer(FINESSE C) 0x0094c23f 0xf1000001

== Detector Buffer(FINESSE A) 0x0094c30d 0x69000001 == Detector Buffer(FINESSE C) 0x0094c30d 0x69000001

0x0094c23f 0xf1000001

0x0094c13a 0x91000001 == Detector Buffer(FINESSE C) 0x0094c13a 0x91000001

```
1, Get dummy data file (data from two CDC FEE boards connected to FINESSE A and C.)
                         login.cc.kek.jp: ~yamadas/rawdata/root output RawCDC rev7133.root
                         2, See contents of the data
                         % cd ${BELLE2 LOCAL DIR}/dag/; svn update
                         % cd ${BELLE2 LOCAL DIR}/daq/rawdata/examples/
                         % basf2 ReadStoreTemplate.py -i ./root output RawCDC rev7133.root | less
[INFO] Steering file: ReadStoreTemplate.py
>>> basf2 Python environment set
>>> Framework object created: fw
                                                                               In this data,
==== DataBlock(RawCDC): Block # 0: Event # 0: node ID 0x00000000: block size 224 bytes
                                                                               Detector buffer contains only 2words(=8bytes)
                                                                               per/FINESSE/event.
==== DataBlock(RawCDC): Block # 1: Event # 1: node ID 0x00000000: block size 224 bytes
                                                                               Note that block # is a number used by DAQ software
                                                                               for handling data and not related with Event #.
==== DataBlock(RawCDC): Block # 2: Event # 2: node ID 0x00000000: block size 224 bytes
```