

RawCOPPER data format

May. 30, 2017 (ff5c3626434)
Satoru Yamada

1, Overview of RawCOPPER format (one data block from a COPPER board)






- RawCOPPER header
 - COPPER header
 - B2link HSLB header (slot A FINNESSE)
 - B2link FEE header(slot A FINNESSE)
 - Data contents(Detector buffer) (slot A FINNESSE)
 - B2link FEE trailer (slot A FINNESSE)
 - B2link HSLB trailer (slot A FINNESSE)
 - B2link HSLB header (slot B FINNESSE)
 - B2link FEE header(slot B FINNESSE)
 - Data contents(Detector buffer) (slot B FINNESSE)
 - B2link FEE trailer (slot B FINNESSE)
 - B2link HSLB trailer (slot B FINNESSE)
 - B2link HSLB header (slot C FINNESSE)
 - B2link FEE header(slot C FINNESSE)
 - Data contents(Detector buffer) (slot C FINNESSE)
 - B2link FEE trailer (slot C FINNESSE)
 - B2link HSLB trailer (slot C FINNESSE)
 - B2link HSLB header (slot D FINNESSE)
 - B2link FEE header(slot D FINNESSE)
 - Data contents(Detector buffer) (slot D FINNESSE)
 - . B2link FEE trailer (slot D FINNESSE)
 - B2link HSLB trailer (slot D FINNESSE)
 - COPPER trailer
- RawCOPPER trailer

HSLB: High speed link board

1-1, Online header/trailer reduction



3

- PreRawCOPPER format
- If you store data by COPPER CPU, then output data will be in Pre(reduction)RawCOPPER format.
- PostRawCOPPER format
- Store the data downstream from readout PC, the output data will be in Post(reduction)RawCOPPERFormat

2-1, "RawCOPPER header/trailer" format in PreRawCOPPER format (ver. 2+0x80)

Use this version number to distinguish * Same as ver.1 Different data format. Ver.0: to 2014. June(including DESY test) ver.1: from June.2014 Ver.2 : from Apr. 2015 (rev.17269) Number of total words 0x7f7f | Format ver.(8bit) | Number of words in this block 2 (8bit) RawCopper Header(20words) exp no. (10bit=1024), run no.(22bit=4194304 including subrun) B2link(FEE event number(32bit) COPPER header (13words) /HSLB) From B2link FEE header 2 (TT-ctime | Trig-type) header From B2link FEE header 2 (TT-utime) and Node ID Detector buffer(slot A) trailer b2link CRC error bit (4) | |truncation mask (truncated or not) / type of data (compressed, calibration, ...) offset to 1st block of user's data(FEE data) offset to 2nd block of user's data(FEE data) Detector buffer(slot B) offset to 3rd block of user's data(FEE data) offset to 4th block of user's data(FEE data) Detector buffer(slot C) Detector buffer(slot D) COPPER trailer(3words?) reserved (maybe used for checksum) RawCopper Trailer(2words) termination word of this block = 0x7fff0006

2-2, "RawCOPPER header" and trailer format in PostRawCOPPER format (ver.0x02)

* Same as ver.1

Same as PreRawCOPPER format

2-3, tentative format of 32bit node ID

Format:

(31-24) Detector ID: 8bit=256: detector ID (9-0) lower bits of COPPER ID: 10bit (1024)

Detector ID:

Detector ID (Defined in rawdata/dataobjects/include/RawCOPPERFormat.h)

- #define SVD ID 0x01000000 // tentative
- #define CDC ID 0x02000000 // tentative
- #define BPID ID 0x03000000 // tentative
- #define EPID ID 0x04000000 // tentative
- #define BECL ID 0x05000000 // tentative
- #define EECL ID 0x06000000 // tentative
- #define BKLM ID 0x07000000 // tentative
- #define EKLM ID 0x08000000 // tentative

- #define TRGDATA ID 0x10000000 // tentative
 - #define CDCTRGDATA ID 0x11000000 // tentative
 - #define TOPTRGDATA ID 0x12000000 // tentative
 - #define ECLTRGDATA ID 0x13000000 // tentative
 - #define KLMTRGDATA ID 0x14000000 // tentative
 - #define GDLTRGDATA ID 0x15000000 // tentative

Full COPPER ID:

```
Except for TRG, full COPPER ID can be reconstructed by "( Detector ID >> 24 ) * 1000 + COPEPR ID(12bit) "
         e.g. NodeID = 0x0600000a -> COPPER ID = cpr6010
            NodeID = 0x0100000a -> COPPER ID = cpr1010
```

- For TRG data, full COPPER ID can be reconstructed by "9000 + COPEPR ID(12bit) "
 - Following a request from SLC part, hostname of TRG coppers will be assigned as follows: (Aug.27, 2015)

CDC TRG: cpr1100* ECL TRG: cpr1200*

GDLTRG: cpr1500*


A label of COPPER ID will be attached on the front of a COPPER board

Node ID = "TTD" = 0x54544420 and "FTSW" are reserved for VME CPU and FTSW now.

3-1, COPPER header and trailer in PreRawCOPPER format (ver. 2 + 0x80)

* Same as ver.1

COPPER header



COPPER Trailer

3-2, COPPER header and trailer in PostRawCOPPER format (ver. 0x02)

* Same as ver.1

No COPEPR header and trailer in Post reduction rawcopper format.

4-1, B2link FEE header/Trailer, B2link HSLB header/Trailer in PreRawCOPPERFormat (ver. 0x02 + 0x80)

* Same as ver.1

From Nakao-san's Belle2link User guide (June 10, 2014): You can download from 18 th B2GM indico page

http://kds.kek.jp/getFile.py/access?contribId=132&sessionId=28&resId=0&materialId=0&confl d=15329




Figure 5: Data format as read out by the COPPER. The header and trailer words labelled with HSL are attached by HSLB, the words with B2L are attached by the belle2link component, and the words with FEE are those written into the belle2link component by the frontend firmware.

NOTICE:

To produce this format, the b2tt core used in

the FEE firmware should be the latest.

Please see Nakao-san's following e-mails:

[b2link_ml:0143] Belle2link version 0.01 -

SVN update

And

[b2link_ml:0144] Re: Belle2link version 0.01 -

SVN update.

4-2, B2link FEE header/Trailer, B2link HSLB header/Trailer in

PostRawCOPPERFormat (ver. 0x02)

FEE: Data #n (32)

B2L: CRC error count(16)

* Modified from ver.1

HSL: 0xFFAA(16) -- B2L header | HSLB-tag(16)

B2L: '0'(1) | TT-ctime(27) | TT-type(4)

B2L: TT-tag(32)

B2L: TT-expran(32)

B2L: '0' | B2L-ctime(27) | reserved(4)

FEE: Data #0 (32)

FEE: Data #1 (32)

FEE:

B2link HSLB header -> 1word

B2link FEE header -> 1word

B2link FEE trailer -> 1word

B2link HSLB trailer -> 0word

Figure 5: Data format as read out by the COPPER. The header and trailer words labelled with HSL are attached by HSLB, the words with B2L are attached by the belle2link component, and the words with FEE are those written into the belle2link component by the frontend firmware.

4-3, Older B2link header/trailer formats

At DESY test in January of 2014

From Nakao-san's B2GM slides:

http://kds.kek.jp/getFile.py/access?contribId=143&sessionId=38&resId=0&materialId=slides&confId=13911

Data format (Final?)

The format used at the telescope test

```
HSL: 0xFFAA(16) --- B2L header | HSLB-tag(16)

B2L: '0'(1) | TT-ctime(27) | TT-type(4)

B2L: TT-tag(32)

B2L: TT-utime(32)

B2L: TT-exprun(32)

B2L: '0' | B2L-ctime(27) | debug-flag(4)

FEE: Data #0 (32)

FEE: Data #1 (32)


FEE: Data #n (32)

B2L: TT-tag(16) | B2L-checksum(16)

HSL: 0xFF55(16) | HSLB checksum(16)
```


tag (event number) and utime to be increased to 32-bit (done),
 HSLB-checksum, B2L-checksum to be added


5-1, RawDataBlock object (to handle Raw data from COPPER board)



Source code:

https://belle2.cc.kek.jp/svn/trunk/software/rawdata/dataobjects/






Example of Data handling on COPPER (as of rev.12453)


RecvSendCOPPER.py

```
DeSerializerCOPPER.cc
                                                                           Serializer.cc
                 Event()
                                                                           Event()
                 // Read 50 evens per an event() loop
                                                                     // Send each
                                                                     raw_dblkarray[i] event
                 For(int I =0; i<
                                                                     separately.
                 NUM_EVT_PER_BASF2LOOP_COPPER(=50);i++){
                                                                     For( int I =0; i<
                 Num_event = 1; num_nodes = 1;
                                                                                                                    Send data
                                                                     raw_dblkarray.getEntries();i++){
                 temp_rawdblk =
                                                                                                                    to a readout
From
                 raw dblkarray.appendNew();
                                                                     sendByWriteV( raw dblkarra
FEE
                                                                                                                    PC
                   temp_rawdblk->SetBuffer(temp_buf,
                                                                     y[ | ])
                 m size word, delete flag, num events,
                 num_nodes)
                 Raw dblkarray.getEntries()=50
                                                                  Raw dblkarray.getEntries()=50
                 Raw dblkarray->GetNumNodes()=1
                                                                  Raw dblkarray->GetNumNodes()=1
                 Raw dblkarray->GetNumEvents()=1
                                                                  Raw dblkarray->GetNumEvents()=1
```

Example of Data handling on a readout PC (as of rev.12453)



Example: # of events and nodes in one RawDataBock(Detector) object and StoreArray:



Example: how to store rawdata in RawDetector object by DAQ program

Full Belle II DAQ

- HLT(High Levele Trigger) receives serialized binary data from readout PCs and stores them in RawDetector Class. This RawDetector objects will be stored in storage.
 - Module : daq/rfarm/event/modules/src/Raw2DsModule.cc

Pocket DAQ

- DATA are stored as a RawDataBlock object
- Example program can be used to convert RawDataBlock/COPPER to RawDetector objects
 - Module : daq/rawdata/modules/src/Convert2RawDet.cc

2-2, Rawdata Unpacker for new and old data formats

Data taken at the DESY beam test(old format) can be read with the latest rawdata package -> by checking data ver. In header.

New RawCOPPER class

- ➤ No change in style of the member functions -> No effect on derived class
- > Does not have a format information in itself
 - > Format class contains format information
 - > RawCOPPERformat.cc -> the latest format
 - RawHeader.cc
 - ➤ RawCOPPERformat v0.cc -> an old format
 - > RawHeader v0.cc
 - Assgin a format class to m_access in CheckVersionSetBuffer()
 - ➤ Use m access to access buffer contents

```
inline int RawCOPPER::GetExpNo(int n)
{
   CheckVersionSetBuffer();
   return m_access->GetExpNo(n);
}
inline int RawCOPPER::GetRunNo(int n)
{
   CheckVersionSetBuffer();
   return m_access->GetRunNo(n);
}
```

Notice:

- RawCOPPER class supports both formats for a while (0.5-1 year after the format becomes stable?).
- In that case, the latest RawCOPPER class cannot be used to read old format
- Of course, you can use old rawdata repository to read old format
- For ver.0 format, use rawdata repository before 11228

Revision History of this document(1)

- Jan.5, 2014 rev. 8376 : Add definition of tentative subsysID format
- Dec. 16, 2013 rev.7974 :
 - Add B2linkFEE header format
 - Add comments about handling StoreArray when unpacking Raw*** data.
- Oct.21, 2013 :rev.7133
 - Add instruction about Rawdata unpacking program
- Oct. 18, 2013 :rev. 7095
 - 1 st draft
- Jun. 23, 2014 : rev. 11234
 - Online (header/trailer) reduction scheme on readout PC is introduced
 - RawHeader format is changed
 - COPPER header/trailer format is changed
 - Nakao-san updated B2LFEE/HSLB header/trailer format
 - See [b2link ml:0144] Re: Belle2link version 0.01 SVN update
- Aug. 23, 2014: rev. 12453
 - Add a description of how RawDataBLock objects are handled by the actual DAQ program.
- Sep. 26, 2014 : rev.13065
 - Add TRG ID definition (0x09000000)
- Oct. 24, 2014 : rev. 13460
 - Add a slide about # of events and nodes in one RawDataBock(Detector) object and StoreArray
 and how to store raw data in RawDetector object by DAQ program (p.16. p.17)
- Jan. 23, 2015: rev. 15030
 - Modify TRGDATA_ID and add ***TRGDATA_ID for trigger from sub-detectors
- Mar.3, 2015 : rev. 15988
 - Use "node ID" instead of subsystem ID.

Revision History of this document(2)

- Apr. 21, 2015 rev. 17269 : Description about a new data format (ver.2)
- Aug. 27, 2015: Hostname of TRG COPPER are revised
- May 30, 2017: NODE_ID for ECLTRG and TOPTRG was swapped. (commit ID: ff5c3626434)

Backup

Modified part of data-format in the change from ver.1 to ver.2

4-2, B2link FEE header/Trailer, B2link HSLB header/Trailer in PostRawCOPPERFormat (ver. 0x01)




Figure 5: Data format as read out by the COPPER. The header and trailer words labelled with HSL are attached by HSLB, the words with B2L are attached by the belle2link component, and the words with FEE are those written into the belle2link component by the frontend firmware.