
Pocket DAQ manual


(2014.8.8) rev.11654 Satoru Yamada

Connection between components



Software components

These are stored under https://belle2.cc.kek.jp/svn/trunk/software/daq



Pocket DAQ w/o slow controller and GUI

O. Before using pocket DAQ

- Setup a PXE boot server for COPPER CPU and install driver for COPPER etc
 - See https://belle2.cc.kek.jp/~twiki/bin/view/Detector/DAQ/PocketDAQ
- Install basf2 on both COPPER CPU and Readout PC
 - <u>See https://belle2.cc.kek.jp/~twiki/bin/view/Software/SoftwareInstallation</u>. "addall" does no install the daq package. You need to do "addpkg daq".
 - Check daq/Sconscripts and check if env['CONTINUE'] = False is commented out.
 - Compile with scons.
- Compile eventbuilder
 - cd \${BELLE2_LOCAL_DIR}/daq/eventbuilder/evb0/; gmake eb0

Install GLEW for running basf2 program

```
After external ver.0.5.0?, basf2 requires GLEW(The OpenGL Extension Wrangler Library).

The following procedure worked fine.

1, install the yum epel repository.

ropc01$ wget http://ftp-srv2.kddilabs.jp/Linux/distributions/fedora/epel/5/i386/epel-release-5-4.noarch.rpm

ropc01$ ssh cpr** -lroot

cpr**# rpm -ivh epel-release-5-4.noarch.rpm

cpr**# ls -lrt /etc/yum.repos.d/

cpr**# exit

2, install glew on COPPER

ropc01$ sudo yum install glew --installroot=/tftpboot/copper/root
```

1, Set parameters(1)

```
[ROPC] % cd ${BELLE2_LOCAL_DIR}/daq/copper/daq_scripts
[ROPC] % emacs run_start.sh
```

Set arguments of start copper.sh

Usage: start_copper.sh <HOSTNAME> <COPPER node> <FINNESSE bit flag: A=1, B=2, C=4, D=8> <Use network shared memory: Yes=1, No=0> <node name(used for Network Shared Memory)> e.g.)

/usr/bin/xterm -fn 7x14 -geometry 102x10+0+642 -e \${BELLE2_LOCAL_DIR}/daq/copper/daq_scripts/start_copper.sh cpr006 1 1 0 0& /usr/bin/xterm -fn 7x14 -geometry 102x10+750+642 -e \${BELLE2_LOCAL_DIR}/daq/copper/daq_scripts/start_copper.sh cpr007 2 3 0 0&

<COPPER node ID> will be attached to RawCOPPER header.

1, Set parameters(2)

NOTICE: Environment variables for basf2

In \${BELLE2_LOCAL_DIR}/daq/copper/daq_scripts/copper.sh, a line, "source ~/.bash_profile", is for setting up basf2 environment. You need to add basf2 setting commands in your .bash_profile (or other script file).

Please see "Setup of Software Tools" at https://belle2.cc.kek.jp/~twiki/bin/view/Software/SoftwareInstallation for details.

```
[ ${BELLE2_LOCAL_DIR}/daq/copper/daq_scripts/copper.sh ]
#
# setup basf2 environment (See "Setup of Software Tools" at
https://belle2.cc.kek.jp/~twiki/bin/view/Computing/SoftwareInstallation
#
source ~/.bash_profile
```

1, Set parameters(3)

Specify whether communication with GUI via shared memory will be used or not:

1, Edit \${BELLE2_LOCAL_DIR}/daq/rawdata/examples/RecvStream1.py

```
If you use GUI, set 'UseShmFlag' = 1
receiver.param('UseShmFlag', 1)
If you don't use GUI set 'UseShmFlag' = 0
receiver.param('UseShmFlag', 0)
```

2, Edit \${BELLE2_LOCAL_DIR}/daq/rawdata/examples/RecvStream1.py

```
If you use GUI, set 'UseShmFlag' = 1
receiver.param('UseShmFlag', 1)
If you don't use GUI set 'UseShmFlag' = 0
receiver.param('UseShmFlag', 0)
```

2, Setting of merger (eb0) (1)

1, compile

```
$ cd ${BELLE2_LOCAL_DIR}/daq/eventbuilder/evb0 $ make
```

2, edit eb0.sh

cd /home/usr/b2daq/eb # change to cd "your local Bellell directroy"/daq/eventbuilder/evb0/eb-xinetd (e.g. cd /home/usr/yamada/basf2/release/daq/eventbuilder/evb0/eb-xinetd)

If you wan to read data from COPPERs with hostnames of cpr010, cpr011 and cpr012, use the following ./eb0 cpr010:33000 cpr011:33000 cpr012:33000

2, Setting of merger (eb0) (2)

3, edit eb0's part in \${BELLE2_LOCAL_DIR}/daq/eventbuilder/evb0/eb-xinted

```
Modify user and server.

service eb0

{
    bind = 127.0.0.1
    socket_type = stream
    type = UNLISTED
    port = 5101
    wait = no
    user = g0cdc # Change to your user name (e.g. yamagata)
    server = /home/usr/b2daq/eb/eb0.sh
    # change to "your Belle2 local directory"/daq/eventbuilder/evb0/eb0.sh

# (e.g. /home/usr/yamagata/basf2/release/daq/eventbuilder/evb0/eb0.sh)
}
```

You need to change the other services in eb-xinetd (eb1tx, eb1rx, eb2tx, eb2rx) in the same way.

2, Setting of merger (eb0) (3)

4, Set link to eb-xinetd

```
# cd /etc/xinetd.d/
# In -s ${BELLE2_LOCAL_DIR}/daq/eventbuilder/evb0/eb-
xinetd eb-xinetd
# /sbin/service xinetd restart
```

The eb0 daemon will be automatically invoked when a basf2 program on a readout PC connects to a port specified in ebxinted.

3, How to start DAQ

[ROPC] % cd \${BELLE2_LOCAL_DIR}/release/daq/copper/daq_scripts [ROPC] % ./run_start.sh

4, How to stop DAQ

- No stop button for now (If you use GUI, you have.)
- You need to specify max # of events or time to stop the run on a basf2 python file.

```
ROPC % cd ${BELLE2_LOCAL_DIR}/daq/rawdata/examples Edit RecvStream1.py
```

You can set following paramters to stop a run. receiver.param('MaxTime', 300.) receiver.param('MaxEventNum', 30.)

-1 means inifinite.

5. Read an output file and extract FEE buffer

- Output file name
 - \${BELLE2_LOCAL_DIR}/daq/copper/daq_scripts/root_output.root
 - You can change the filename by editing \${BELLE2_LOCAL_DIR}/daq/rawdata/examples/RecvStream1.py
- Read the root file
 - Please see dataformat and unpacker manuals.
 - Data format
 - https://belle2.cc.kek.jp/svn/trunk/software/daq/copper/doc/SetupPocketDAQ_4p1_Raw COPPERDataHandling.pdf
 - Unpacker
 - https://belle2.cc.kek.jp/svn/trunk/software/daq/copper/doc/SetupPocketDAQ_4p5_Raw DataUnpackerPacker.pdf

Revision history

- 2013/10/16
 - Change explanation about option of run_start.sh script
- 2014/7/15 (rev. 11619)
 - To specify FINESSE bit field, you need to use decimal value (not hexadecimal).
 - Add an argument of RecvSenDCOPPER.py
 - Update eb0 related issues.
 - About eb0.sh
 - About eb0-xinted
- 2014/7/16(rev. 11654)
 - eb0.sh and eb-xinetd need more modification
- 2014/7/19
 - Add a comment about "addpkg daq"
- 2014/8/8
 - Modify usage of start_run.sh. (nodename is only used for NSM)

Appendix

A-1, yum setting for COPPER

Sometime yum fails when you try to install software to COPPER directory; readoutPC # yum install **** --installroot=/tftpboot/copper/root

In that case, you need to change yum setting for COPPER host:

- 1, Copy yum repository files. readoutPC # cp /etc/yum.repo.d/* /tftpboot/copper/root/etc/yum.repo.d/*
- 2, Set a proxy server in /tftpboot/copper/root/etc/yum.conf same as readoutPC:/etc/yum.conf

A-2, Socket parameters on COPPER, readout PC and etc..

- 1, Enable "Jumbo Frame" in network switch's setting, if necessary (e.g. HP1810-24G)

 If the switch does not support Jumbo frame, setting MTU a large value may cause a problem.
- 2, Change socket parameters
- #/sbin/ifconfig eth? mtu 9000

```
# /sbin/sysctl -w net.ipv4.tcp_rmem="8388608 8388608 8388608"
# /sbin/sysctl -w net.ipv4.tcp_wmem="8388608 8388608 8388608"
# /sbin/sysctl -w net.ipv4.tcp_mem="9000000 9000000 9000000"
# /sbin/sysctl -w net.core.rmem_default=8388608
# /sbin/sysctl -w net.core.wmem_default=8388608
# /sbin/sysctl -w net.core.rmem_max=8388608
# /sbin/sysctl -w net.core.wmem_max=8388608
```

3, edit /etc/rc.local(or /tftpboot/copper/root/etc/rc.local) if you want the parameters automatically set at boot time

A-3, Use ntpd to adjust time on COPPER

cpr**** \$ chkconfig –list | grep ntp

A-4, syslogd

Install syslogd:

readoutPC # yum install syslogd --installroot=/tftpboot/copper/root

Add entry to /tftpboot/copper/root/etc/syslogd.conf if needed:

Log anything (except mail) of level info or higher.

Don't log private authentication messages!

*.info;mail.none;news.none;authpriv.none;cron.none /var/log/messages

B-1, Test bench at Tsukuba B3

